



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MRPL47 Antibody [#C14083]

Catalog Number: C14083

Concentration: 1mg/ml

Swiss-Prot No.: Q9HD33

Other Names: 39S ribosomal protein L47, mitochondrial; nasopharyngeal carcinoma metastasis-related 1; NCM1; RM47

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MRPL47.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MRPL47 antibody detects endogenous levels of total MRPL47 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Muzny D.M., Nature 440:1194-1198(2006).

Mural R.J., Submitted (SEP-2005) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team; Genome Res. 14:2121-2127(2004).

