



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ATPG Antibody [#C14596]

Catalog Number: C14596

Concentration: 1mg/ml

Swiss-Prot No.: P36542

Other Names: ATP synthase gamma chain, mitochondrial precursor; ATP5C1; ATPG; EC 3.6.3.14

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ATPG.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ATPG antibody detects endogenous levels of total ATPG protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Matsuda C., J. Biol. Chem. 268:24950-24958(1993).

Deloukas P., Nature 429:375-381(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

