



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ATP5G2 Antibody [#C14598]

Catalog Number: C14598

Concentration: 1mg/ml

Swiss-Prot No.: Q06055

Other Names: ATP synthase lipid-binding protein, mitochondrial; ATP synthase proteolipid P2; ATPase protein 9; ATPase subunit c

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ATP5G2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ATP5G2 antibody detects endogenous levels of total ATP5G2 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:40000

References:

Dyer M.R., Biochem. J. 293:51-64(1993).

Higuti T., Biochim. Biophys. Acta 1173:87-90(1993).

Farrell L.B., Biochem. Biophys. Res. Commun. 144:1257-1264(1987).



Peptide -

+

Immunohistochemistry analysis of paraffin-embedded human pancreas tissue using ATP5G2 antibody.