



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ATPAF2 Antibody [#C14600]

Catalog Number: C14600

Concentration: 1mg/ml

Swiss-Prot No.: Q8N5M1

Other Names: ATP synthase mitochondrial F1 complex assembly factor 2; ATP12; ATPF2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ATPAF2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ATPAF2 antibody detects endogenous levels of total ATPAF2 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Wan D., Proc. Natl. Acad. Sci. U.S.A. 101:15724-15729(2004).

Ota T., Nat. Genet. 36:40-45(2004).

Zody M.C., Nature 440:1045-1049(2006).

