



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ATP5S Antibody [#C14607]

Catalog Number: C14607

Concentration: 1mg/ml

Swiss-Prot No.: Q99766

Other Names: ATP synthase subunit s, mitochondrial [Precursor]; ATP synthase-coupling factor B; Mitochondrial ATP synthase regulatory component factor B; ATP5S; ATPW

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ATP5S.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ATP5S antibody detects endogenous levels of total ATP5S protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Belogradov G.I., J. Biol. Chem. 277:6097-6103(2002).

Ota T., Nat. Genet. 36:40-45(2004).

Mural R.J., Submitted (SEP-2005) to the EMBL/GenBank/DDBJ databases.

