



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ABCA6 Antibody [#C14615]

Catalog Number: C14615

Concentration: 1mg/ml

Swiss-Prot No.: Q8N139

Other Names: ABC transporter ABCA6; ABCA6;
ATP-binding cassette A6; ATP-binding cassette sub-family
A member 6; ATP-binding cassette, sub-family A (ABC1),
member 6; EST155051; FLJ43498

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from C-terminal of human
ABCA6.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: ABCA6 antibody detects endogenous levels
of total ABCA6 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Kaminski W.E., Biochem. Biophys. Res. Commun.
285:1295-1301(2001).

Arnould I., GeneScreen 1:157-164(2001).

Ota T., Nat. Genet. 36:40-45(2004).



For Research Use Only