



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ABCA8 Antibody [#C14617]

Catalog Number: C14617

Concentration: 1mg/ml

Swiss-Prot No.: O94911

Other Names: ATP-binding cassette sub-family A member 8

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ABCA8.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ABCA8 antibody detects endogenous levels of total ABCA8 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Nagase T., DNA Res. 5:355-364(1998).

Ohara O., Submitted (JAN-2004) to the EMBL/GenBank/DDBJ databases.

Tsuruoka S., Biochem. Biophys. Res. Commun. 298:41-45(2002).

