



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ABCD1 Antibody [#C14624]

Catalog Number: C14624

Concentration: 1mg/ml

Swiss-Prot No.: P33897

Other Names: ABC42; adrenoleukodystrophy protein; ALD; ALDP; AMN; ATP-binding cassette, sub-family D (ALD), 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ABCD1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ABCD1 antibody detects endogenous levels of total ABCD1 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Mosser J., Nature 361:726-730(1993).

Ross M.T., Nature 434:325-337(2005).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

