



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ABCD4 Antibody [#C14627]

Catalog Number: C14627

Concentration: 1mg/ml

Swiss-Prot No.: O14678

Other Names: ATP-binding cassette sub-family D member 4; Peroxisomal membrane protein 69; PMP69; Peroxisomal membrane protein 1-like; PXMP1-L; P70R

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ABCD4.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ABCD4 antibody detects endogenous levels of total ABCD4 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)
Mouse (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:10000

References:

Farmer A., Submitted (MAY-2003).

Korn B., Submitted (JUN-2004).

Nakamura Y., J. Hum. Genet. 47:285-310(2002).

