



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

HTR2C Antibody [#G008]

Catalog Number: G008

Concentration: 1mg/ml

Swiss-Prot No.: P28335

Other Names: 5-hydroxytryptamine 2C receptor; 5H2C; 5HT-1C; 5HT1C; 5HT2C; HTR1C; HTR2C; Serotonin receptor

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human HTR2C.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: HTR2C antibody detects endogenous levels of total HTR2C protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Saltzman A.G., Biochem. Biophys. Res. Commun. 181:1469-1478(1991).

Stam N.J., Eur. J. Pharmacol. 269:339-348(1994).

Xie E., Genomics 35:551-561(1996).

