



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

HTR5A Antibody [#G015]

Catalog Number: G015

Concentration: 1mg/ml

Swiss-Prot No.: P47898

Other Names: 5-hydroxytryptamine receptor 5A; 5-HT-5A;
5-HT-5; Serotonin receptor 5A

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human HTR5A.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: HTR5A antibody detects endogenous levels of total HTR5A protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:


Rees S., FEBS Lett. 355:242-246(1994).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Sjoebloom T., Science 314:268-274(2006).



Western blot analysis of extracts from Jurkat cells, using HTR5A antibody.



Immunofluorescence analysis of LOVO cells, using HTR5A antibody.