



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## JMY Antibody [#C10003]

**Catalog Number:** C10003

**Concentration:** 1mg/ml

**Swiss-Prot No.:** A1L4P5

**Other Names:** JMY protein, JMY

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human JMY.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** JMY antibody detects endogenous levels of total JMY protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:1000

**References:**

Strausberg R.L., Proc. Natl. Acad. Sci. U.S.A. 99:16899-16903(2002).

NIH MGC Project, Submitted (JAN-2007) to the EMBL/GenBank/DDBJ databases.

