



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## TTF2 Antibody [#C10022]

**Catalog Number:** C10022

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O00358

**Other Names:** Forkhead box protein E1, Thyroid transcription factor 2, TTF-2, Forkhead-related protein FKHL15, FKHL15, TTF2, TTF2

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human TTF2.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** TTF2 antibody detects endogenous levels of total TTF2 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Chadwick B.P., Genomics 41:390-396(1997).

Macchia P.E., Biochimie 81:433-440(1999).

Baris I., J. Clin. Endocrinol. Metab. 91:4183-4187(2006).

