



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## CDKA1 Antibody [#C10037]

**Catalog Number:** C10037

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O14519

**Other Names:** Cyclin-dependent kinase 2-associated protein 1, CDK2-associated protein 1, Putative oral cancer suppressor, Deleted in oral cancer 1, DOC-1, CDK2AP1, CDKAP1, DOC1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human CDKA1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** CDKA1 antibody detects endogenous levels of total CDKA1 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:1000

**References:**

Daigo Y., Genes Chromosomes Cancer 20:204-207(1997).

Tsuji T., J. Biol. Chem. 273:6704-6709(1998).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

