



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

DLX3 Antibody [#C10136]

Catalog Number: C10136

Concentration: 1mg/ml

Swiss-Prot No.: O60479

Other Names: Homeobox protein DLX-3, DLX3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human DLX3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: DLX3 antibody detects endogenous levels of total DLX3 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100


ELISA: 1:10000

References:

Price J.A., Hum. Mol. Genet. 7:563-569(1998).


Sumiyama K., Proc. Natl. Acad. Sci. U.S.A. 99:780-785(2002).

Dong J., Am. J. Med. Genet. A 133:138-141(2005).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human heart tissue using DLX3 antibody.



Western blot analysis of extracts from 293 cells, using DLX3 antibody.