



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MEKKK 4 Antibody [#C10242]

Catalog Number: C10242

Concentration: 1mg/ml

Swiss-Prot No.: O95819

Other Names: EC 2.7.11.1, HPK/GCK-like kinase HGK, M4K4, MAP4K4, MAPK/ERK kinase kinase 4, MEK kinase kinase 4, MEKKK 4, Nck interacting kinase, hepatocyte progenitor kinase-like/germinal center kinase-like kinase, mitogen-activated protein kinase kinase kinase 4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human MEKKK 4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MEKKK 4 antibody detects endogenous levels of total MEKKK 4 protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:5000

References:


Yao Z., J. Biol. Chem. 274:2118-2125(1999).

Wright J.H., Mol. Cell. Biol. 23:2068-2082(2003).

Ishikawa K., DNA Res. 5:169-176(1998).



Western blot analysis of extracts from COLO cells and HUVEC cells, using MEKKK 4 antibody.



Immunofluorescence analysis of A549 cells, using MEKKK 4 antibody.