



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

DGKQ Antibody [#C10713]

Catalog Number: C10713

Concentration: 1mg/ml

Swiss-Prot No.: P52824

Other Names: DAG kinase theta, DAGK, DAGK4, DAGK7, DGK-theta, diacylglycerol kinase theta, diglyceride kinase theta

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human DGKQ.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: DGKQ antibody detects endogenous levels of total DGKQ protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Pilz A., Genomics 26:599-601(1995).

