



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

UBXD5 Antibody [#C11186]

Catalog Number: C11186

Concentration: 1mg/ml

Swiss-Prot No.: Q5T124

Other Names: UBX domain-containing protein 5, Socius, Colorectal tumor-associated antigen COA-1, SOC, PP2243

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human UBXD5.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: UBXD5 antibody detects endogenous levels of total UBXD5 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100


ELISA: 1:20000

References:

Maccalli C., Cancer Res. 63:6735-6743(2003).

Ota T., Nat. Genet. 36:40-45(2004).


Wan D., Proc. Natl. Acad. Sci. U.S.A. 101:15724-15729(2004).



Peptide -

+

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using UBXD5 antibody.



Peptide -

+

Western blot analysis of extracts from HepG2 cells, using UBXD5 antibody.