



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TTC23 Antibody [#C11200]

Catalog Number: C11200

Concentration: 1mg/ml

Swiss-Prot No.: Q5W5X9

Other Names: TPR repeat protein 23, Tetratricopeptide repeat protein 23, Cervical cancer proto-oncogene 8 protein, HCC-8, HCC8

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TTC23.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TTC23 antibody detects endogenous levels of total TTC23 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Kim J.W., Submitted (AUG-2001) to the EMBL/GenBank/DDBJ databases.

Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

