



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PKNOX2 Antibody [#C11510]

Catalog Number: C11510

Concentration: 1mg/ml

Swiss-Prot No.: Q96KN3

Other Names: Homeobox protein PKNOX2, PBX/knotted homeobox 2, Homeobox protein PREP-2, PKNOX2, PREP2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human PKNOX2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PKNOX2 antibody detects endogenous levels of total PKNOX2 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Imoto I., Biochem. Biophys. Res. Commun. 287:270-276(2001).

Fognani C., Nucleic Acids Res. 30:2043-2051(2002).

