



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ATF6B Antibody [#C11585]

Catalog Number: C11585

Concentration: 1mg/ml

Swiss-Prot No.: Q99941

Other Names: Cyclic AMP-dependent transcription factor ATF-6 beta, Activating transcription factor 6 beta, ATF6-beta, cAMP-responsive element-binding protein-like 1, cAMP response element-binding protein-related protein, Creb-rp, Protein G13, Processed cyclic AMP-dependent transcription factor ATF-6 beta, CREBL1, G13

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human ATF6B.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ATF6B antibody detects endogenous levels of total ATF6B protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Min J., Genomics 30:149-156(1995).

Khanna A., Biochem. J. 319:81-89(1996).

Xie T., Genome Res. 13:2621-2636(2003).

