



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## NRIP2 Antibody [#C11598]

**Catalog Number:** C11598

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9BQI9

**Other Names:** Nuclear receptor-interacting protein 2, NRIP2

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from N-terminal of human NRIP2.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** NRIP2 antibody detects endogenous levels of total NRIP2 protein.


**Reactivity:** Human (Identities = 100%, Positives = 100%);  
Mouse (Identities = 81%, Positives = 100%)

**Applications:** IHC: 1:50~1:100 ELISA: 1:20000

**References:**

Wiemann S., Genome Res. 11:422-435(2001).

The MGC Project Team, Genome Res. 14:2121-2127(2004).



Peptide                      -                      +

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using NRIP2 antibody.