



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ADCK4 Antibody [#C11484]

Catalog Number: C11484

Concentration: 1mg/ml

Swiss-Prot No.: Q96D53

Other Names: EC 2.7.11.- ; ADCK4; Uncharacterized aarF domain-containing protein kinase 4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ADCK4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: ADCK4 antibody detects endogenous levels of total ADCK4 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:

The MGC Project Team, Genome Res. 14:2121-2127(2004).
Ota T., Nat. Genet. 36:40-45(2004).



Western blot analysis of extracts from HepG2 cells, using ADCK4 antibody.



Immunofluorescence analysis of A549 cells, using ADCK4 antibody.