



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MASTL Antibody [#C11491]

Catalog Number: C11491

Concentration: 1mg/ml

Swiss-Prot No.: Q96GX5

Other Names: Microtubule-associated serine/threonine-protein kinase-like, EC 2.7.11.1, THC2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human MASTL.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MASTL antibody detects endogenous levels of total MASTL protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Ota T., Nat. Genet. 36:40-45(2004).

Deloukas P., Nature 429:375-381(2004).

Beausoleil S.A., Nat. Biotechnol. 24:1285-1292(2006).

