



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ALS2CR11 Antibody [#C14459]

Catalog Number: C14459

Concentration: 1mg/ml

Swiss-Prot No.: Q53TS8

Other Names: Amyotrophic lateral sclerosis 2
chromosomal region candidate gene 11 protein

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from C-terminal of human
ALS2CR11.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: ALS2CR11 antibody detects endogenous
levels of total ALS2CR11 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Ota T., Nat. Genet. 36:40-45(2004).

Hillier L.W., Nature 434:724-731(2005).

Sjoebloom T., Science 314:268-274(2006).

