



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ALS2CR4 Antibody [#C14462]

Catalog Number: C14462

Concentration: 1mg/ml

Swiss-Prot No.: Q96Q45

Other Names: ALS2S4; amyotrophic lateral sclerosis 2 chromosomal region candidate gene 4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ALS2CR4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ALS2CR4 antibody detects endogenous levels of total ALS2CR4 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Hadano S., Nat. Genet. 29:166-173(2001).

Ota T., Nat. Genet. 36:40-45(2004). Hillier L.W., Nature 434:724-731(2005).

