



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

KAL1 Antibody [#C14510]

Catalog Number: C14510

Concentration: 1mg/ml

Swiss-Prot No.: P23352

Other Names: adhesion molecule-like X-linked; ADMLX; anosmin-1; HHA; KAL; KALIG-1; KALIG1; Kallmann syndrome 1; KALM; KMS

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human KAL1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KAL1 antibody detects endogenous levels of total KAL1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Legouis R., Cell 67:423-435(1991).

del Castillo I., Nat. Genet. 2:305-310(1992).

Franco B., Nature 353:529-536(1991).

