



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## APOA5 Antibody [#C14525]

**Catalog Number:** C14525

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q6Q788

**Other Names:** Apolipoprotein A-V; Apo-AV; ApoA-V;  
Apolipoprotein A5; Regeneration-associated protein 3

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human APOA5.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** APOA5 antibody detects endogenous levels of total APOA5 protein.

**Reactivity:** Human


**Applications:** WB: 1:500~1:3000 ELISA: 1:1000

**References:**

van Der Vliet H.N., J. Biol. Chem. 276:44512-44520(2001) [PubMed: 11577099].

Beckstead J.A., Biochemistry 42:9416-9423(2003) [PubMed: 12899628].

Fullerton S.M., Hum. Genet. 115:36-56(2004) [PubMed: 15108119].



**For Research Use Only**