



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

APOF Antibody [#C14533]

Catalog Number: C14533

Concentration: 1mg/ml

Swiss-Prot No.: Q13790

Other Names: Apolipoprotein F; Apo-F; Lipid transfer inhibitor protein; LTIP

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human APOF.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: APOF antibody detects endogenous levels of total APOF protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Day J.R., Biochem. Biophys. Res. Commun. 203:1146-1151(1994).

Wang X., J. Biol. Chem. 274:1814-1820(1999).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

