



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ARG2 Antibody [#C14560]

Catalog Number: C14560

Concentration: 1mg/ml

Swiss-Prot No.: P78540

Other Names: Arginase-2 mitochondrial; Arginase II; Non-hepatic arginase; Kidney-type arginase

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human ARG2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ARG2 antibody detects endogenous levels of total ARG2 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Gotoh T., FEBS Lett. 395:119-122(1996).

Vockley J.G., Genomics 38:118-123(1996).

Morris S.M. Jr., Gene 193:157-161(1997).

