



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ARSD Antibody [#C14566]

Catalog Number: C14566

Concentration: 1mg/ml

Swiss-Prot No.: P51689

Other Names: Arylsulfatase D; EC 3.1.6.-; ASD

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ARSD.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ARSD antibody detects endogenous levels of total ARSD protein.


Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IHC: 1:50~1:100 IF: 1:100~1:500
ELISA: 1:20000

References:


Franco B., Cell 81:15-25(1995).

Urbitsch P., DNA Cell Biol. 19:765-773(2000).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using ARSD antibody.



Peptide - +

Immunofluorescence analysis of MCF-7 cells, using ARSD antibody.