



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

LW-1 Antibody [#C11810]

Catalog Number: C11810

Concentration: 1mg/ml

Swiss-Prot No.: Q9UBD0

Other Names: Heat shock transcription factor family X-linked 1, LW-1, HSFX1, hCG_27154

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human LW-1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: LW-1 antibody detects endogenous levels of total LW-1 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Strausberg R.L., Proc. Natl. Acad. Sci. U.S.A. 99:16899-16903(2002).

Wang L., Submitted (APR-1999) to the EMBL/GenBank/DDBJ databases.

Mural R.J., Submitted (JUL-2005) to the EMBL/GenBank/DDBJ databases.

