



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## DBF4 Antibody [#C11821]

**Catalog Number:** C11821

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9UBU7

**Other Names:** Protein DBF4 homolog A, Activator of S phase kinase, DBF4-type zinc finger-containing protein 1, Chiffon homolog A, DBF4, ASK, DBF4A, ZDBF1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from N-terminal of human

DBF4.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** DBF4 antibody detects endogenous levels of total DBF4 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Kumagai H., Mol. Cell. Biol. 19:5083-5095(1999).

Jiang W., EMBO J. 18:5703-5713(1999).

Ota T., Nat. Genet. 36:40-45(2004).

