



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

NXF1 Antibody [#C11823]

Catalog Number: C11823

Concentration: 1mg/ml

Swiss-Prot No.: Q9UBU9

Other Names: MEX67, TAP, Tip-associated protein, mRNA export factor TAP, nuclear RNA export factor 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human NXF1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: NXF1 antibody detects endogenous levels of total NXF1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Braun I.C., EMBO J. 18:1953-1965(1999).

Kang Y., Genes Dev. 13:1126-1139(1999).

Bear J., Mol. Cell. Biol. 19:6306-6317(1999).

