



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## Mlx Antibody [#C11841]

**Catalog Number:** C11841

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9UH92

**Other Names:** BigMax protein, MAD7, MXD7, Max-like bHLHZip protein, Max-like protein X, TCFL4, protein Mlx, transcription factor-like 4, transcription factor-like protein 4

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human Mlx.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** Mlx antibody detects endogenous levels of total Mlx protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Billin A.N., J. Biol. Chem. 274:36344-36350(1999).

Meroni G., Oncogene 19:3266-3277(2000).

Ota T., Nat. Genet. 36:40-45(2004).

