



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TRPS1 Antibody [#C11843]

Catalog Number: C11843

Concentration: 1mg/ml

Swiss-Prot No.: Q9UHF7

Other Names: Zinc finger transcription factor Trps1, Zinc finger protein GC79, Tricho-rhino-phalangeal syndrome type I protein, TRPS1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human TRPS1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TRPS1 antibody detects endogenous levels of total TRPS1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:1000

References:

Momeni P., Nat. Genet. 24:71-74(2000).

Chang G.T.G., J. Natl. Cancer Inst. 92:1414-1421(2000).

Olsen J.V., Cell 127:635-648(2006).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using TRPS1 antibody.