



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PAR4 (Cleaved-Gly48) Antibody [#L0377]

Catalog Number: L0377

Concentration: 1mg/ml

Swiss-Prot No.: Q96R10

Other Names: Proteinase-activated receptor 4
[Precursor], PAR-4, Thrombin receptor-like 3, Coagulation
factor II receptor-like 3, F2RL3, PAR4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from human PAR4.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: PAR4 (Cleaved-Gly48) antibody detects
endogenous levels of fragment of activated PAR4
resulting from cleavage adjacent to Gly48.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Kahn M.L., J. Biol. Chem. 273:23290-23296(1998).

Kahn M.L., Nature 394:690-694(1998).

Xu W.-F., Proc. Natl. Acad. Sci. U.S.A.
95:6642-6646(1998).

