



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

SUMO2/3 (Cleaved-Gly93) Antibody [#L0379]

Catalog Number: L0379

Concentration: 1mg/ml

Swiss-Prot No.: P61956

Other Names: HSMT3; MGC117191; SMT3 (suppressor of mif two 3, yeast) homolog 2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human SUMO2/3.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SUMO2/3 (Cleaved-Gly93) antibody detects endogenous levels of fragment of activated SUMO2/3 resulting from cleavage adjacent to Gly93.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Mannen H., Biochem. Biophys. Res. Commun. 222:178-180(1996).

Lapenta V., Genomics 40:362-367(1997).

Kamitani T., J. Biol. Chem. 273:11349-11353(1998).

