



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

THRB (AP2, Cleaved-Arg327) Antibody [#L0383]

Catalog Number: L0383

Concentration: 1mg/ml

Swiss-Prot No.: P00734

Other Names: Coagulation factor II; EC 3.4.21.5; F2;
THRB

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human THRB.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: THRB (AP2, Cleaved-Arg327) antibody detects endogenous levels of fragment of activated THRB resulting from cleavage adjacent to Arg327.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Degen S.J.F., Biochemistry 26:6165-6177(1987).

Degen S.J.F., Biochemistry 22:2087-2097(1983).

Walz D.A., Proc. Natl. Acad. Sci. U.S.A. 74:1969-1972(1977).

