



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## FOXD3 Antibody [#C11862]

**Catalog Number:** C11862

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9UJU5

**Other Names:** Forkhead box protein D3, HNF3/FH transcription factor genesis, FOXD3

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human FOXD3.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** FOXD3 antibody detects endogenous levels of total FOXD3 protein.

**Reactivity:** Human, Mouse


**Applications:** WB: 1:500~1:3000 IHC: 1:50~1:100  
ELISA: 1:5000

**References:**

Gregory S.G., Nature 441:315-321(2006).


Guo Y., Proc. Natl. Acad. Sci. U.S.A. 99:3663-3667(2002).

Ramsey H.C., Submitted (OCT-1999) to the EMBL/GenBank/DDBJ databases.



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human heart tissue using FOXD3 antibody.



Western blot analysis of extracts from HUVEC cells and 293 cells, using FOXD3 antibody.

**For Research Use Only**