



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MAFF Antibody [#C11887]

Catalog Number: C11887

Concentration: 1mg/ml

Swiss-Prot No.: Q9ULX9

Other Names: Transcription factor MafF, V-maf musculoaponeurotic fibrosarcoma oncogene homolog F, U-Maf, MAFF

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MAFF.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MAFF antibody detects endogenous levels of total MAFF protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:10000

References:

Ota T., Nat. Genet. 36:40-45(2004).

Moran J.A., Biochem. J. 361:371-377(2002).

Dunham I., Nature 402:489-495(1999).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using MAFF antibody.