



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

MMP17 (Cleaved-Gln129) Antibody [#L0325]

Catalog Number: L0325

Concentration: 1mg/ml

Swiss-Prot No.: Q9ULZ9

Other Names: EC 3.4.24.-; matrix metalloproteinase-17; MMP-17; MT-MMP 4; MT4-MMP; MT4MMP

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human MMP17.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: MMP17 (Cleaved-Gln129) antibody detects endogenous levels of fragment of activated MMP17 resulting from cleavage adjacent to Gln129.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Kajita M., FEBS Lett. 457:353-356(1999).

Seiki M., Submitted (JAN-2002) to the EMBL/GenBank/DDBJ databases.

Puente X.S., Cancer Res. 56:944-949(1996).

