



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## MMP27 (Cleaved-Tyr99) Antibody [#L0334]

**Catalog Number:** L0334

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9H306

**Other Names:** Matrix metalloproteinase-27 [Precursor];  
EC 3.4.24.-; MMP-27; MMP27; UNQ2503/PRO5992

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from human MMP27.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** MMP27 (Cleaved-Tyr99) antibody detects endogenous levels of fragment of activated MMP27 resulting from cleavage adjacent to Tyr99.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Benoit de Coignac A., Submitted (OCT-1999) to the EMBL/GenBank/DDBJ databases.

Clark H.F., Genome Res. 13:2265-2270(2003).

Bar-Or A., Brain 126:2738-2749(2003).

