



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## NOTCH2 (Cleaved-Ala1734) Antibody [#L0352]

**Catalog Number:** L0352

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q04721

**Other Names:** Neurogenic locus notch homolog protein 2 [Precursor], Notch 2, hN2, Notch 2 extracellular truncation, Notch 2 intracellular domain, NOTCH2

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from human NOTCH2.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** NOTCH2 (Cleaved-Ala1734) antibody detects endogenous levels of fragment of activated NOTCH2 resulting from cleavage adjacent to Ala1734.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Stifani S., Nat. Genet. 2:119-127(1992).

Blaumueller C.M., Cell 90:281-291(1997).

Gray G.E., Am. J. Pathol. 154:785-794(1999).

