



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## ZNF337 Antibody [#C11951]

**Catalog Number:** C11951

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9Y3M9

**Other Names:** Zinc finger protein 337, ZN337, ZNF337

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human ZNF337.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ZNF337 antibody detects endogenous levels of total ZNF337 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Deloukas P., Nature 414:865-871(2001).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

The German cDNA consortium; Submitted (MAY-1999) to the EMBL/GenBank/DDBJ databases.

