



Web: www.anbio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TBX22 Antibody [#C11954]

Catalog Number: C11954

Concentration: 1mg/ml

Swiss-Prot No.: Q9Y458

Other Names: T-box transcription factor TBX22, T-box protein 22, TBX22, TBOX22

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human TBX22.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TBX22 antibody detects endogenous levels of total TBX22 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Braybrook C., Nat. Genet. 29:179-183(2001).

Ross M.T., Nature 434:325-337(2005).

Laugier-Anfossi F., Gene 255:289-296(2000).

