



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

SPIN1 Antibody [#C11983]

Catalog Number: C11983

Concentration: 1mg/ml

Swiss-Prot No.: Q9Y657

Other Names: OCR, Ovarian cancer-related protein, Spindlin-1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human SPIN1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SPIN1 antibody detects endogenous levels of total SPIN1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Zhang H.L., Chin. Sci. Bull. 45:909-914(2000).

Gao Y., Biochem. Biophys. Res. Commun. 335:343-350(2005).

Humphray S.J., Nature 429:369-374(2004).

