



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

14-3-3 β Antibody [#C12001]

Catalog Number: C12001

Concentration: 1mg/ml

Swiss-Prot No.: P31946

Other Names: 14-3-3 protein beta/alpha, 1433B, 143B, KCIP-1, Protein 1054, Protein kinase C inhibitor protein-1, YWHAB

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human 14-3-3 β .

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: 14-3-3 β antibody detects endogenous levels of total 14-3-3 β protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Leffers H., J. Mol. Biol. 231:982-998(1993).

Deloukas P., Nature 414:865-871(2001).

Gevaert K., Nat. Biotechnol. 21:566-569(2003).

Chi A., J. Proteome Res. 5:3135-3144(2006).

