



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

14-3-3 ϵ Antibody [#C12002]

Catalog Number: C12002

Concentration: 1mg/ml

Swiss-Prot No.: P62258

Other Names: 14-3-3E, 143E, KCIP-1, Mitochondrial import stimulation factor L subunit, Protein kinase C inhibitor protein-1, YWHAE

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human 14-3-3 ϵ .

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: 14-3-3 ϵ antibody detects endogenous levels of total 14-3-3 ϵ protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Conklin D.S., Proc. Natl. Acad. Sci. U.S.A. 92:7892-7896(1995).

Jin D.-Y., Nature 382:308-308(1996).

Chong S.S., J. Biol. Chem. 271:735-741(1996).

