



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

5-HT-1A Antibody [#C12008]

Catalog Number: C12008

Concentration: 1mg/ml

Swiss-Prot No.: P08908

Other Names: 5-HT1A, 5-hydroxytryptamine 1A receptor, 5HT1A, G-21, HTR1A, Serotonin receptor 1A

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human 5-HT-1A.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: 5-HT-1A antibody detects endogenous levels of total 5-HT-1A protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Kobilka B.K., Nature 329:75-79(1987).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Parks C.L., Nucleic Acids Res. 19:7155-7160(1991).

